

# **HYDROGEN PEROXIDE**

# **Determination by Ceric Sulfate Method**

**Technical Data Sheet** 

## DETERMINATION of HYDROGEN PEROXIDE (H<sub>2</sub>O<sub>2</sub>): CERIC SULFATE METHOD

### **SCOPE**

This method is suitable for the determination of hydrogen peroxide in the presence of other peroxygens. The range covered is 5 - 200 g/kg of hydrogen peroxide.

## **PRINCIPLE**

The sample is diluted in dilute sulfuric acid cooled with ice. The hydrogen peroxide is titrated with ceric sulfate solution using ferroin as an indicator.

Caution: Safety glasses must be worn when using this method.

#### REAGENTS

All reagents should be analytical reagent grade unless otherwise stated. Distilled water, or water of equivalent purity, should be used in this method.

Ceric sulfate solution (0.1 N).

Sulfuric acid solution (1:19) v/v.

Carefully and slowly add 50 mL of concentrated sulfuric acid (d 1.84 g/mL) to 950 mL of water in a 2-liter beaker with continuous stirring and cooling.

Caution: Safety goggles must be worn when handling concentrated sulfuric acid.

Ferroin indicator solution

Dissolve 0.1740 g of ferrous sulfate in 25 mL of water. Add 0.3712 g 1,10 - phenanthroline hydrate to the solution and stir to dissolve. Alternatively a commercially available ferroin indicator solution may be used.

## **APPARATUS**

Weighing bottles, glass or polyethylene, 1-5 mL capacity.

## **PROCEDURE**

- 1. Weigh (to ± 0.001 g) about 0.4 g of sample into a weighing bottle. Let the mass of the sample taken be Wg.
- 2. Measure 150 mL of sulfuric acid (1:19) into a 500-mL conical flask and add sufficient crushed ice to maintain a temperature of less than 10°C during the titration. When the temperature of the solution is below 10°C, add about three drops of ferroin indicator solution and titrate drop-wise with ceric sulfate solution (0.1N) until the indicator changes to a blue color.
- 3. Add the weighing bottle containing the sample to the cold solution and swirl to mix. Titrate rapidly with ceric sulfate solution (0.1N) to the same blue color.
- 4. Let the titration obtained be A mL.





# **HYDROGEN PEROXIDE**

# **Determination by Ceric Sulfate Method**

**Technical Data Sheet** 

### **CALCULATION**

Hydrogen peroxide content (as  $H_2O_2$ )= A x N x 17.

 $\frac{1 \times 17.01 \times 1000 \text{ g/kg}}{1000 \times W} = \frac{A \times N \times 17.01 \text{ g/kg}}{W}$ 

Where N is the Normality of the ceric sulfate solution.

### **FIRST AID**

In case of product splashing into the eyes and face, treat eyes first.

Eye contact: Flush eyes immediately with water for at least 15 minutes. Call a physician.

**Skin contact:** Immediately flush skin with water while removing contaminated clothing and shoes. Call a physician if irritation persists.

**Inhalation:** Remove the victim from the contaminated area to fresh air. Call a physician in case of respiratory symptoms.

**Ingestion:** Consult with a physician immediately in all cases. DO NOT induce vomiting. If victim is conscious, rinse mouth and give fresh water.

### **DANGER**

Hydrogen peroxide solutions are strong oxidizers and corrosive to the eyes, mucous membranes and skin. Consult the SDS for the appropriate Personal Protective Equipment to wear when handling hydrogen peroxide. In case of contact with the eyes, skin or clothing, flush with large amounts of water for 15 minutes. In case of ingestion, sit upright, drink large quantities of water to dilute the stomach contents and seek immediate medical attention. Product in contact with combustible materials may cause fires.

Before using, read Safety Data Sheet (SDS) for this chemical. Solvay Chemicals, Inc. 24-hour Emergency Phone Number – 800-424-9300 (CHEMTREC®)

To our actual knowledge, the information contained herein is accurate as of the date of this document. However, neither Solvay Chemicals, Inc., nor any of its affiliates, makes any warranty, express or implied, or accepts any liability in connection with this information or its use. This information is for use by technically skilled persons at their own discretion and risk and does not relate to the use of this product in combination with any other substance or any other process. This is not a license under any patent or other proprietary right. The user alone must finally determine suitability of any information or material for any contemplated use in compliance with applicable law, the manner of use and whether any patents are infringed. This information gives typical properties only and is not to be used for specification purposes. Solvay Chemicals, Inc. reserves the right to make additions, deletions or modifications to the information at any time without prior notification.

Trademarks: Trademarks and/or other Solvay Chemicals, Inc. products referenced herein are either trademarks or registered trademarks of Solvay Chemicals, Inc. or its affiliates, unless otherwise indicated.

Before using, read the Safety Data Sheet (SDS) for the chemical.

